

Research Report

Treatment of occlusal caries with CURODONT™ REPAIR and Fluoride

Authors

University of Greifswald
Department of Preventive and Pediatric Dentistry
17475 Greifswald, Germany

Mohammad Alkilzy, Ahmad Tarabaih, Christian Splieth



Treatment of occlusal caries with CURODONT™ REPAIR and Fluoride

Safety, application and clinical effect of CURODONT™ REPAIR in children with early occlusal caries

Summary

Occlusal surfaces of erupting permanent molars are highly prone to caries. The self-assembling peptide (P11-4) has been proven to enhance biomimetic mineralization of early carious lesions.

The aim of this study was to evaluate safety, clinical applicability and effect of using P11-4 (CURODONT™ REPAIR) in non-invasive treatment of early occlusal lesions.

Method

70 patients (28 females, mean age 10.03 years \pm 2.7, dft 2.8 \pm 3.1, DMFT 1.3 \pm 2.5) with early occlusal lesions (ICDAS-II:1-3) on first or second permanent molars at eruption were allocated in this randomized, controlled, single blinded post-marketing study either to test

(CURODONT™ REPAIR and Duraphat®) or control (Duraphat®). Safety and applicability was evaluated using dentist's / patient's questionnaires about adverse events, difficulties of application and satisfaction with the procedure. Lesions were assessed at baseline and recalls after 3 and 6 months regarding clinical status (ICDAS-II), caries activity and DIAGNOdent. At every recall fluoride varnish was applied and patients received oral health instructions.

Results

Preliminary data showed good patient acceptance for CURODONT™ REPAIR. Investigators considered the application as much easier as a composite filling or even a fissure sealant. In all cases, no adverse events or allergic reactions have been observed after application.

Study size

Control group:
35 treated with Duraphat® (Fluoride 22,000 ppm)

Test group:
35 treated with CURODONT™ REPAIR and Duraphat®

Main Selection criteria

- 1) Occlusal caries (ICDAS 1 or 2) on a freshly erupted molar
- 2) No need for imminent interventional treatment
- 3) Informed consent of parent / guardian and patient

Study design

Randomised, assessor blinded, gold standard controlled clinical trial

Diagnostic

DIAGNOdent; ICDAS-II activity status (visual)

Treatment protocol

- 1) Professional dental hygiene treatment
- 2) 2% Sodium-Hypochloride (20 s),
to remove the pellicle
- 3) Etching Gel (35% Phosphoric Acid; 20 s)
to remove amorphous mineral from the pores
to the subsurface lesion
- 4) Neutralising with water (20 s)
- 5) Test group:
Application of CURODONT™ REPAIR, (3–5 min)
- 6) Test- and control grup:
Application of Duraphat®

Follow-up procedure

Follow-up D90 (Additional application of Duraphat®)

Follow-up D180 (Last visit)

Result

| Criteria | Duraphat® | CURODONT™ REPAIR & Duraphat® |
|--|--|--|
| DIAGNODent change D ₀ to D ₉₀ | -4.4 | -13.0 (p=0.06) |
| to D ₁₈₀ | -1.1 | -18.6 (p=0.005) |
| Caries activity | D ₀ : 97 % D ₉₀ : 80 % D ₁₈₀ : 65 % | D ₀ : 100 % D ₉₀ : 46 % (OR=4.7; p=0.02) D ₁₈₀ : 20 % (OR=7.6; p<0.001) |

Conclusion

Duraphat® was effective in stopping the progression of caries. However, the additional treatment of CURODONT™ REPAIR (followed by Duraphat®), led to a significant decrease in lesion progression and its activity.

Literature

Alkilzy et al. (2014): 61st ORCA, July 2–5, 2014, Greifswald, Germany; Caries Research, (2014), Abstract #61
Alkilzy et al. (2015): 62nd ORCA, July 1-4, 2015, Brussels, Belgium; Caries Research, (2015), Abstract #30
www.credentis.com; www.curodont.com

EURODONT™ REPAIR

Natural filling of initial carious lesions
and similar defects



EURODONT™ PROTECT

Regenerating tooth gel for intensive prophylaxis

For further information,
please visit our Website:

www.curodont.com

